

PRODUCT CATALOGUE

health products

About us

Pharm Line Group was founded in 1997 as a representative of a British pharmaceutical company PharmLine Ltd.

The Company specializes on hospital products, such as: syringes, infusion and transfusion systems, vide assortment of wound dressings and bandages from traditional to high-technological, plasters, ovulation and pregnancy tests for early diagnostics. The products are manufactured in Great Britain according to the international quality standards ISO 9002, ISO 13485 µ CE.



Starting from 2005 Pharm Line Group is the official representative of the Japanese Okamoto Industries in Europe, the manufacturer of the thinnest and the most durable condoms in the word.











- The best price
- The wound dressing doesn't stick to the wound due to the special a traumatic coverage
- Hypoallergenic polyacrylic adhesive securely fixes plaster but tears off without pain







- perforated skin color polymer base
- water-resistant
- breathable
- bactericidal (rivanol)
- 3 size in kit



SENSITIVE





- The best offer from the point of view of price and quality ratio
- The wound dressing doesn't stick to the wound due to the special a traumatic coverage
- Hypoallergenic polyacrylic adhesive securely fixes plaster but tears off without pain

20 pcs.

ELASTIC



- 100% cotton fabric base
- breathable
- elastic
- bactericidal (rivanol)
- 3 size in kit



TRAVEL KIT



- perforated skin color polymer base
- water-resistant
- breathable
- bactericidal (rivanol)
- 2 size in kit



SENSITIVE



- non woven base
- super-breathable
- tears off without pain
- bactericidal (rivanol)
- 1 size in kit





- The best offer from the point of view of price and quality ratio
- The wound dressing doesn't stick to the wound due to the special a traumatic coverage
- Hypoallergenic polyacrylic adhesive securely fixes plaster but tears off without pain

INVISIBLE



- transparent polymer base
- water-resistant
- breathable
- bactericidal (rivanol)
- 1 size in kit



- transparent polyurethane film – the "second skin" technology
- 100% waterproof protects the wound while bathing
- 2 size in kit



KIDS



- polymer base with drawings
- water-resistant
- breathable
- bactericidal (rivanol)
- 3 colors in kit





- The best quality plasters in boxes with fist opening control system
- The wound dressing doesn't stick to the wound due to the special a traumatic coverage
- Hypoallergenic polyacrylic adhesive securely fixes plaster but tears off without pain

ELASTIC



- 100% cotton fabric base
- breathable
- elastic
- bactericidal (silver ions)



UNIVERSAL

- perforated skin color polymer base
- water-resistant
- breathable
- bactericidal (silver ions)



SENSITIVE





С ионами серебра



- The best quality plasters in boxes with fist opening control system
- The wound dressing doesn't stick to the wound due to the special a traumatic coverage
- Hypoallergenic polyacrylic adhesive securely fixes plaster but tears off without pain

COMFORT FLEX



- super-elastic polymer base
- water-resistant
- breathable
- bactericidal (silver ions)



AQUA



- transparent polyurethane film the "second skin" technology
- 100% waterproof protects the wound while bathing



KIDS



- polymer base with drawings
- water-resistant
- breathable
- bactericidal (silver ions)



Blister plasters

HYDROCOLLOID



- Hydrocolloid plasters create optimal conditions for fast regeneration of skin
- Lass on the skin for more than 24 hours
- Protects from friction and pressure



Corn plasters

DISSOLVE DRY CORNS



- Effective way to remove dry corns
- Traditional means in a convenient form •







Pepper plaster



- Effective remedy for rheumatic and muscular pain, arthritis, radiculitis, lumbago, colds
- Provides deep heat and increased blood circulation



stops the pain

Fixation tapes

Secure bandages, compresses, catheters...



- fabric tape made of 100% cotton
- elastic material
- zinc oxide adhesive based on natural rubber
- high adhesion



- non-woven base made of microporous cellulose
- breathable
- hypoallergenic polyacrylic adhesive
- can be easily torn by hand
- tears off without pain
- recommended for children and people with sensitive skin



MASTER UNI

Fixation tapes

TENER'S[®]

sizes 5x500 sm 4x500 sm

3x500 sm

2x500 sm

1x500 sm

Secure bandages, compresses, catheters...



- fabric tape made of 100% cotton
- elastic material
- zinc oxide adhesive based on natural rubber – high adhesion
- cellular coating of the adhesive retained breathability of the fabric



FENER'S

T-FIX

3 X 500 cm

- breathable
- hypoallergenic polyacrylic adhesive

SENSTIVE

- can be easily torn by hand
- tears off without pain
- recommended for children and people with sensitive skin



POLYMER

- innovative product
- breathable

TEMERIS

3 X 500 cm

FENER'S

- water-resistant
- hypoallergenic polyacrylic adhesive
- can be easily torn by hand



Pregnancy tests



Principle of the test is based on detecting the "pregnancy hormone" HCG (human chorionic gonadotropin) in the urine. In contact with the hormone chemical reaction takes place and particles of a litmus in the test change color. In the early stages of pregnancy the concentration of HCG (the pregnancy hormone) in your body increases very fast and after 1-2 weeks after conception is 25-100 mME/ml. Conception occurs at the time of ovulation – i.e. approximately in the middle of the menstrual cycle. Thus by the time the anticipated start of the menstruation pregnancy is already about 2 weeks old and the concentration of the hormone is high enough for detection.



Pregnancy tests

FEMITEST

type/ quantity in box

stripe-test/1 pcs.

stripe-test/2 pcs.

jet /1 pcs.

- accuracy of laboratory analyses (accuracy more than 99.5%)
- definition of pregnancy at home from the 1st day of delay
- the sensitivity of 20 mMe/ml manufactured in the UK



The most popular and simple type of aptitude test is a pregnancy test strip. Principle of the test is based on detecting the "pregnancy hormone" HCG (human chorionic gonadotropin) in the urine. In contact with the hormone chemical reaction takes place and particles of a litmus in the test change color. This test is not difficult – just immerse the test vertically in the container with urine to the level indicated by the arrows for 10 seconds, and then place the strip on a dry and horizontal surface and after 3-5 minutes evaluate the results.

The jet tests are more complex: by a fibrous core consisting of a plurality of tubules, the fluid quickly and uniformly rises to a test part where a reagent colors upon contact with HCG. To use the inkjet test is very easy, you don't even need to collect urine in a separate cup – just substitute test, tip down, under a stream of urine, place the test on a horizontal surface and evaluate the results in 3-5 minutes.

name

express

double control

expert



Pregnancy tests

- accuracy of laboratory analyses (accuracy more than 99.5%)
- very early diagnosis the definition of pregnancy 4 days before the expected date of menstruation
- sensitivity 10 mME/ml
- manufactured in the UK STRIPE-TEST



nametype/ quantity in boxultrastripe-test/ 1 pcs.ultra expertjet /1 pcs.





The most popular and simple type of aptitude test is a pregnancy test strip. Principle of the test is based on detecting the "pregnancy hormone" HCG (human chorionic gonadotropin) in the urine. In contact with the hormone chemical reaction takes place and particles of a litmus in the test change color. This test is not difficult – just immerse the test vertically in the container with urine to the level indicated by the arrows for 10 seconds, and then place the strip on a dry and horizontal surface and after 3-5 minutes evaluate the results.

The jet tests are more complex: by a fibrous core consisting of a plurality of tubules, the fluid quickly and uniformly rises to a test part where a reagent colors upon contact with HCG. To use the inkjet test is very easy, you don't even need to collect urine in a separate cup – just substitute test, tip down, under a stream of urine, place the test on a horizontal surface and evaluate the results in 3-5 minutes.

50

TOOC

TOOC

T

T C C



FEMiTEST ultra

Ovulation tests

determine 2 the most favorable day for conception

STRIPE-TESTS

- accuracy of laboratory analyses (accuracy more than 99.5%)
- manufactured in the UK

JET TESTS O serverun **FEMIPLAN FEMIPLAN** expert НА ОВУЛЯЦИЮ O annua труйных тестов C (D starman 5 stripe-tests 5 jet tests

The most popular and simple type of ovulation test is the strip-test. The principle of operation of the test is based on the detection of LH (luteinizing hormone). The highest concentration of LH in the body is indicated 48 hours prior to the release of an egg from the follicle (ovulation). In contact with the hormone chemical reaction takes place and particles of a litmus in the test change color. The test should be done 5 days in a row, to detect the presence of LH and its sufficient concentration. When two strips manifest clearly - the concentrations of LH reached a peak, and within 2 days there will come ovulation. These 2 days are the most favorable days for conception.

The jet tests are more complex: on the fibrous core consisting of a plurality of tubules, the fluid quickly and uniformly rises to a test part where a reagent colors upon contact with LH (luteinizing hormone). To use the inkjet test is very easy, you don't even need to collect urine in a separate cup – just substitute test, tip down, under a stream of urine, place the test on a horizontal surface and evaluate the results in 3-5 minutes.



FEMIPLAN

name	type/ quantity in box
FemiPlan	stripe-test/ 5 pcs.
FemiPlan expert	jet /5 pcs.

(no need to collect urine)

Condoms



- The finest and most durable latex condom in the world (listed in the Guinness Book of Records)
- Manufactured in Japan

shape

№ 1 sales in Japan, USA, Germany, France, Norway (ten times the winner of the annual World Best Condoms prize) •



- skin color
- skin color

classic shape

Condoms



- Thin and durable condoms with abundant lubricant
- Manufactured in UK from premium quality latex







3-component syringes



1 ml with double insulin tuberculin scale

needle diameter according to ISO:

type	diameter	length
	0,45 mm	12 mm
	0,63 mm	32 mm
	0,7 mm	40 mm
	0,8 mm	40 mm

volume	needle	quantity in box
50 ml	0,8*40	5 pcs.
20 ml	0,8*40	5 pcs.
10 ml	0,8*40	5 pcs.
5 ml	0,7*40	10 pcs.
3 ml	0,63*32	10 pcs.
2 ml	0,63*32	10 pcs.
1 ml ins.	0,45*12	10 pcs.

- rubber gasket from synthetic rubber doesn't contain latex, which eliminates completely the risk of allergic reactions
- special silicone lubricant provides a more uniform motion of the piston and provides smooth infusion of medicines without jerks



5-10 syringes in carton box

particularly important for a slow jet injections, accurate dosing of drugs in anesthesia and intensive care





- the cylinder of the syringe is completely transparent it makes possible to control the air bubbles (to exclude embolism of small blood vessels) and to evaluate the proper form of the drug
- a clear contrast indelible scale permits a very precise dosage of the medication
- the retaining ring prevents accidental slipping of the piston from the cylinder
- needles from superior quality steel with triple sharpened tip, silicone treated (unlike regular bevel, edge-triangular sharpening minimizes pain during the injection, silicone reduces skin trauma during punching)
- the needle is slipped onto the cylinder



Terms of cooperation

Our company is distinguished by an individual approach and affordable pricing for each of our client:

- ✓ System of bonuses and discounts from the base price list for distributors and wholesalers
- The possibility of a delay in payment
- Shipment of goods to all regions. We will promptly deliver the goods to the recommended vehicle of your organization
- Professional marketing and trade marketing support
- Commercial and trading equipment
- ✓ The possibility of Contract manufacturing under Your private label

Our team of managers will provide You with expert help in selecting the optimal assortment and the conditions most favorable for You.

Will be glad to answer all your questions :

Contacts:

PharmLine Ltd., Great Britain Cornwall Buildings, 45-51 Newhall Str., office 330, Birmingham, B3 3QR, Great Britain +44 121 2134701

e-mail: <u>info@pharmline.eu</u> web: <u>www. pharmline.eu</u>

